

Anti-Covid-19 Spike antibody



Model	STJ193235
Host	Humanized B Cell
Reactivity	Covid-19
Applications	ECL, ELISA, Immunogold, WB
Immunogen	The recombinant SARS-CoV-2 (2019-nCoV) Spike Protein S
Dilution range	ELISA 1:5000-1:100000WB 1:1000-1:10000
Specificity	This antibody detects SARS-CoV-2 (2019-nCoV) Spike Protein S1 domain
Purification	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Clone ID	S03FC
Note	FOR RESEARCH USE ONLY (RUO).
Clonality	Monoclonal
Conjugation	Unconjugated
Isotype	IgM
Formulation	PBS
Concentration	0.1mg/mL
Storage Instruction	Store at -20°C, and avoid repeat freeze-thaw cycles.